

The Straight Forearm Plate





Straight Forearm Plate
is indicated for fixation of Ulna/Radius
midshaft fractures.

Locking Technology

The Locking plate provides excellent stability on the osteoporotic bone, severely comminuted fractures etc.



Low Profile

The holes and screw heads are designed to enable low profile screw/bone/plate interface. It minimizes the soft tissue irritation risk after the surgery.



Limited Contact

The plate and bone interface minimized to reduce constriction of the blood supply to the periosteum.



Compression Holes

Allow compression and/or distraction.



Screws

The available screws for the Locking Forearm Plates are 3.5mm locking and non-locking cortical screws and 4.0mm cancellous screws.

Ordering Information

| Part No. | Description |
|------------|--|
| 10-1010-04 | 3.5mm Straight Locking Plate, 4 Hole, with 1 Compression Hole |
| 10-1010-06 | 3.5mm Straight Locking Plate, 6 Hole, with 2 Compression Hole |
| 10-1010-08 | 3.5mm Straight Locking Plate, 8 Hole, with 2 Compression Hole |
| 10-1010-10 | 3.5mm Straight Locking Plate, 10 Hole, with 2 Compression Hole |
| 10-1010-12 | 3.5mm Straight Locking Plate, 12 Hole, with 4 Compression Hole |
| 10-1010-14 | 3.5mm Straight Locking Plate, 14 Hole, with 4 Compression Hole |

SCREWS

| | |
|------------|---|
| 30-1271-18 | 2.7mm x 18.0mm Cortical Screw, Self Tapping |
| 30-1271-20 | 2.7mm x 20.0mm Cortical Screw, Self Tapping |
| 30-1271-22 | 2.7mm x 22.0mm Cortical Screw, Self Tapping |
| 30-1271-24 | 2.7mm x 24.0mm Cortical Screw, Self Tapping |
| 30-1271-26 | 2.7mm x 26.0mm Cortical Screw, Self Tapping |
| 30-1271-28 | 2.7mm x 27.5mm Cortical Screw, Self Tapping |
| 30-1271-30 | 2.7mm x 30.0mm Cortical Screw, Self Tapping |
| 30-1271-32 | 2.7mm x 32.0mm Cortical Screw, Self Tapping |
| 30-1271-34 | 2.7mm x 34.0mm Cortical Screw, Self Tapping |
| 30-1271-36 | 2.7mm x 36.0mm Cortical Screw, Self Tapping |
| 30-1271-38 | 2.7mm x 38.0mm Cortical Screw, Self Tapping |
| 30-1271-40 | 2.7mm x 40.0mm Cortical Screw, Self Tapping |
| 30-1271-45 | 2.7mm x 45.0mm Cortical Screw, Self Tapping |
| 30-1271-50 | 2.7mm x 50.0mm Cortical Screw, Self Tapping |
| 30-1271-55 | 2.7mm x 55.0mm Cortical Screw, Self Tapping |
| 30-1271-60 | 2.7mm x 60.0mm Cortical Screw, Self Tapping |
| 30-1271-65 | 2.7mm x 65.0mm Cortical Screw, Self Tapping |

| | |
|------------|---|
| 30-2271-12 | 2.7mm x 12.0mm Locking Cortical Screw, Self Tapping |
| 30-2271-14 | 2.7mm x 14.0mm Locking Cortical Screw, Self Tapping |
| 30-2271-16 | 2.7mm x 16.0mm Locking Cortical Screw, Self Tapping |
| 30-2271-18 | 2.7mm x 18.0mm Locking Cortical Screw, Self Tapping |
| 30-2271-20 | 2.7mm x 20.0mm Locking Cortical Screw, Self Tapping |
| 30-2271-22 | 2.7mm x 22.0mm Locking Cortical Screw, Self Tapping |

| | |
|------------|---|
| 30-1351-08 | 3.5mm x 8.0mm Cortical Screw, Self Tapping |
| 30-1351-10 | 3.5mm x 10.0mm Cortical Screw, Self Tapping |
| 30-1351-12 | 3.5mm x 12.0mm Cortical Screw, Self Tapping |
| 30-1351-14 | 3.5mm x 14.0mm Cortical Screw, Self Tapping |
| 30-1351-16 | 3.5mm x 16.0mm Cortical Screw, Self Tapping |
| 30-1351-18 | 3.5mm x 18.0mm Cortical Screw, Self Tapping |
| 30-1351-20 | 3.5mm x 20.0mm Cortical Screw, Self Tapping |
| 30-1351-22 | 3.5mm x 22.0mm Cortical Screw, Self Tapping |
| 30-1351-24 | 3.5mm x 24.0mm Cortical Screw, Self Tapping |
| 30-1351-25 | 3.5mm x 25.0mm Cortical Screw, Self Tapping |
| 30-1351-26 | 3.5mm x 26.0mm Cortical Screw, Self Tapping |
| 30-1351-28 | 3.5mm x 27.5mm Cortical Screw, Self Tapping |
| 30-1351-30 | 3.5mm x 30.0mm Cortical Screw, Self Tapping |
| 30-1351-32 | 3.5mm x 32.0mm Cortical Screw, Self Tapping |
| 30-1351-34 | 3.5mm x 34.0mm Cortical Screw, Self Tapping |

| Part No. | Description |
|------------|---|
| 30-1351-36 | 3.5mm x 36.0mm Cortical Screw, Self Tapping |
| 30-1351-38 | 3.5mm x 38.0mm Cortical Screw, Self Tapping |
| 30-1351-40 | 3.5mm x 40.0mm Cortical Screw, Self Tapping |
| 30-1351-45 | 3.5mm x 45.0mm Cortical Screw, Self Tapping |
| 30-1351-50 | 3.5mm x 50.0mm Cortical Screw, Self Tapping |
| 30-1351-55 | 3.5mm x 55.0mm Cortical Screw, Self Tapping |
| 30-1351-60 | 3.5mm x 60.0mm Cortical Screw, Self Tapping |
| 30-1351-65 | 3.5mm x 65.0mm Cortical Screw, Self Tapping |
| 30-2351-08 | 3.5mm x 8.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-10 | 3.5mm x 10.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-12 | 3.5mm x 12.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-14 | 3.5mm x 14.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-16 | 3.5mm x 16.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-18 | 3.5mm x 18.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-20 | 3.5mm x 20.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-22 | 3.5mm x 22.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-24 | 3.5mm x 24.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-26 | 3.5mm x 26.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-28 | 3.5mm x 28.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-30 | 3.5mm x 30.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-32 | 3.5mm x 32.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-34 | 3.5mm x 34.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-36 | 3.5mm x 36.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-38 | 3.5mm x 38.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-40 | 3.5mm x 40.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-45 | 3.5mm x 45.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-50 | 3.5mm x 50.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-55 | 3.5mm x 55.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-60 | 3.5mm x 60.0mm Locking Cortical Screw, Self Tapping |
| 30-2351-65 | 3.5mm x 65.0mm Locking Cortical Screw, Self Tapping |
| 30-1401-12 | 4.0mm x 12.0mm Cancellous Screw, Self Tapping |
| 30-1401-14 | 4.0mm x 14.0mm Cancellous Screw, Self Tapping |
| 30-1401-16 | 4.0mm x 16.0mm Cancellous Screw, Self Tapping |
| 30-1401-18 | 4.0mm x 18.0mm Cancellous Screw, Self Tapping |
| 30-1401-20 | 4.0mm x 20.0mm Cancellous Screw, Self Tapping |
| 30-1401-22 | 4.0mm x 22.0mm Cancellous Screw, Self Tapping |
| 30-1401-24 | 4.0mm x 24.0mm Cancellous Screw, Self Tapping |
| 30-1401-26 | 4.0mm x 26.0mm Cancellous Screw, Self Tapping |
| 30-1401-28 | 4.0mm x 28.0mm Cancellous Screw, Self Tapping |
| 30-1401-30 | 4.0mm x 30.0mm Cancellous Screw, Self Tapping |
| 30-1401-32 | 4.0mm x 32.0mm Cancellous Screw, Self Tapping |
| 30-1401-34 | 4.0mm x 34.0mm Cancellous Screw, Self Tapping |
| 30-1401-35 | 4.0mm x 35.0mm Cancellous Screw, Self Tapping |
| 30-1401-36 | 4.0mm x 36.0mm Cancellous Screw, Self Tapping |
| 30-1401-38 | 4.0mm x 38.0mm Cancellous Screw, Self Tapping |
| 30-1401-40 | 4.0mm x 40.0mm Cancellous Screw, Self Tapping |
| 30-1401-45 | 4.0mm x 45.0mm Cancellous Screw, Self Tapping |
| 30-1401-50 | 4.0mm x 50.0mm Cancellous Screw, Self Tapping |
| 30-1401-55 | 4.0mm x 55.0mm Cancellous Screw, Self Tapping |
| 30-1401-60 | 4.0mm x 60.0mm Cancellous Screw, Self Tapping |

The Straight Forearm Plates

manufactured by



www.responseortho.com

Response Ortho LLC
725 River Road, Suite 32-254
Edgewater NJ 07024, USA
Phone: +1 (201) 203 5773
Fax: +1 (201) 608 7169

CE 1015